

Final Project Review

Stride
April 25, 2018



Group Members



Advisor
Prof. Tessier



Richie Hartnett
CSE



Joe Menzie
CSE



Jarred Penney
EE



Jack Higgins
EE

What is the Problem?




- Parkinson's Disease (PD) makes walking challenging
- Physical therapy and other treatments are expensive
- Limited inexpensive methods of monitoring exercises outside of clinical environment

What is Stride?

- Low cost array of wearable sensors that collects body movement information, designed for those with Parkinson's Disease
- Provide real-time feedback and track long term performance progress
- Used in home as well as in clinical environment

Neuro Training


- Vibrations on knee to indicate stride length error
 - Patient-specific stride length threshold
 - Two levels of vibrations depending on severity of error (both benign)
- Auditory Cueing
 - Rhythms played via app to stimulate proper cadence
 - Recommended by a Neurological PT



System Requirements


- Accurately collect movement data to appropriately monitor an individual's:
 - Stride length
 - Cadence
 - Heel-to-toe motion
 - Freezing
- Provide real-time feedback to correct stride length during exercise (less than 100ms)
- Store and display data Android application to track long term patient progress
- Lightweight product that is easy for patient to put on
 - Sensor systems < 1 pound each
 - Waist clip (Raspberry Pi + power supply) < 1 pound
- Sensor systems and Raspberry Pi will have battery life of greater than 2 hours

FPR Block Diagram



Hardware Components (Sensor Systems)

Hardware Component	Shoe (x2)	Knee Sleeve (x2)
PCB	Yes	Yes
Processor	Atmega328p	Atmega328p
Wifi	No	No
Bluetooth	HC-05	HC-05
Battery	Lithium Polymer	Lithium Polymer
Feedback Circuit	No	Yes
IMU	BNO055	BNO055



Hardware Components (Waist Clip Pouch)


Waist Clip Pouch

- Raspberry Pi
 - Processor
 - Bluetooth
 - Wifi
- Lithium Ion Battery




Data Analysis-Ankle

- Outputs: stride length and cadence
 - How it works:
 - Determine step by examining sensor orientation
 - Capture smallest angle in step = θ
 - Stride Length (theoretical) = $L \cdot \cos(\theta)$, where L is length of leg
 - Stride Length (after experimentation) = $L \cdot (K) \cdot e^{(-M \cdot \theta)}$
- (L = leg length; K, M = calibration constants;
 θ = angular extent of swing phase)
- Cadence = steps/time



Data Analysis-Ankle (Cont'd)



Trendline of measured vs. calculated stride length


- Stride length equation and feedback work best when calibrated to the specific user
- Can set feedback thresholds based on severity of disease; therapist's discretion
- Use progressive values or standardized values

Data Analysis-Shoe

- Outputs: heel-to-toe weight transfer and freezes
- How it works:
 - Utilize orientation data to determine range of angles in steps
 - Greater range of angle = better heel-to-toe weight transfer
 - Freeze recorded if angle does not change above 5° threshold





Visual System Communication Overview



System Communication: Bluetooth Connections


- HC-05
- Start and Stop signals
- Synchronous intake



System Communication: MQTT Connections



- Publish/Subscribe Messaging
- App commands
- Amazon IoT Topic

```
dataTopic    Dec 8, 2017 10:41:44 AM -0500    Export  Hide
{
  "sessionId": "1.49",
  "magnetometer": "0.66",
  "accelerometer x": "-3.18",
  "user": "rhartnett",
  "accelerometer z": "8.61",
  "magnetometer x": "0.26",
  "magnetometer y": "1.03",
  "gyroscope x": "-0.42",
  "gyroscope y": "2.03",
  "gyroscope z": "-7.84",
  "accelerometer y": "-2.02"
}
```



System Communication: Other Wifi Connections

- Application
- Rules are used to relay information



Application Backend

- Used AWS Mobile Hub to integrate backend with Android Application
- Created multiple NoSQL database tables to store application data
- Database consists of 6 tables
 - User Information
 - Patient List for Therapist
 - Patient Workout List
 - Overall Workout Performance
 - Session Workout Performance
 - Session Error Check



Application



- Android application designed in Android Studio
- Functionality:
 - Patient/Therapist Login
 - Therapist can program workouts for specific patient
 - Therapist can view patients session performance
 - Patient can view session performance
 - Patient can view and perform programmed workouts
 - Patient workout auditory queuing

Proposed FDR Deliverables

- Fully functional knee and shoe sensor with data analysis
- Breadboard → PCB
- Real-time feedback incorporated into data analysis
- Fully functional and polished application
 - Be able to start workout from application

Goals of Demo

- Show basic data analysis using knee sensor system and view data received by the Raspberry Pi over BlueTooth
- Show data movement through all systems
 - BNO055 → Atmega328 → HC-05 → Raspberry Pi → AWS backend → Android Application
- Show simultaneous sensor data transfer
- Show feedback system functionality
- Show application design and functionality
 - login/create account
 - Therapist side
 - Patient side

Thank You

Questions?